



This is a picture of the
video EEG cart
that will be in your room.




10 uV Sens 0.53 Hz LF 70 Hz HF IA/DB OFF Pat Ref

Elapsed 00:00:00 Remaining 10:00:00

EEG Storage Drive C




Sony Network Camera... http://192.168.1.110



Navigation controls: arrows, WIDE, TELE, CLOSE

The camera is controlled by using a mouse. The video EEG tech will show whoever is staying with you how to use it.

On the next page you can see how the screen will look; an EEG screen with your picture on it! You can also see the screens for controlling the camera.





These two cameras are on the cart, one is for daytime pictures and the other is for night time pictures (infrared). The daytime pictures are in color and the nighttime pictures are in black and white.